INSTRUCTION FOR USE



NOCOSPRAY 2 REMOTE CONTROL

The remote control enable to switch NOCOSPRAY® 2 device on and off, from outside a room.

The documentation applies to Nocospray 2.0 remote control in its original manufacturer's condition which has not been altered.

Review this manual before using the remote control. It was written to ensure you master all functionalities.

Characteristics

Item reference: 3010.611

Wireless technogoly

Frequency: 2,4 GHz

Battery: CR2032

Voltage: 3V

Reach: 20 metres when no obstacle





Battery replacement



- 1 Remove the screw.
- 2 Open the housing.
- 3 Remove the old battery.
- 4 Dispose of the battery in suitable collection points.
- 5 Install the new battery by pressing it. Check the polarity (« + » face up)
- 5 Close the housing.
- 6 Insert the screw and tighten it.

OXY'PHARM

829 rue Marcel Paul 94500 Champigny-sur-Marne commercial@oxypharm.net T: +33145187870





Use of the remote control

Pairing with Nocospray (« blank » remote control)

All remote controls are delivered « blank » and must be paired to a Nocospray® device to be working.

Before first use, follow the instructions below:



Nocospray needs to be turned off but plugged in



Press remote control button (+/- 5 seconds) remaining very close to the Nocospray device that needs to be paired until the Led flashes



While maintaining pressure on the remote control button, switch the Nocospray that needs to be paired on (both actions at the same time)



Remote control Led starts flashing quickly and the paired Nocospray bleeps several times to signal that pairing was successful

Pairing to another Nocospray (« Re-pairing »)



The remote control can be paired to only one Nocospray® at a time.

However, a remote control already paired can be « re-paired » to another Nocospray®

« Re-pairing » method :



Check that all Nocospray close to the remote control are well switched off including the one to be re-paired



Plug in the Nocospray to be re-paired



Press long on remote control button (+/- 30 secondes) remaining very close to the Nocospray device that needs to be repaired until the Led flashes











During these 30 seconds, the Led light flashes and stops alternately.
This phase is repeated 3 times.

Once the Led light on the remote control flashes regularly and without stopping, switch the Nocospray that needs to be re-paired on, while maintaining pressure on the remote control button (both actions at the same time)

Remote control Led starts flashing quickly and the paired Nocospray bleeps several times to signal that pairing was successful

Recycling

IMPORTANT INFORMATION ABOUT THE ELIMINATION COMPATIBLE WITH THE ENVIRONMENT

This device is subject to separate collection for electrical and electronic equipment in accordance with Directive 2002/96 / EC (applicable in European Union countries and other European countries with separate collection systems).

This symbol that appears on the device indicates that at the end of its life cycle, the device will need to be treated separately from household waste.

The device will be assigned to a center of selective collection for electrical and electronic appliances, or reported to the dealer on the purchase of a new device.



The user is responsible for the delivery of the device at the end of life to the collection authorities or to a corporation duly empowered under penalty of sanctions provided for by the legislation on the elimination of waste. The waste supplies must be managed in compliance with the local regulations relating to the management of waste.

Proper differential collection, and the subsequent recycling, processing and environmentally compatible disposable of waste equipment avoids unnecessary damage to the environment and possible related health risks and also promotes recycling of the materials used in the appliance