



Toulouse, December 7th 2017

STUDY 17-2217

Determination of sporicidal activity for aerial surface disinfection processes
According to the method described
in the standard *NF T 72-281* (November 2014)

Human Health

Additional Conditions : *Clostridium difficile* (spore)

Client

OXY'PHARM
917 rue Marcel Paul
94500 CHAMPIGNY SUR MARNE

Test laboratory

FONDEREPHAR
Faculté des Sciences Pharmaceutiques
35 Chemin des Maraîchers
31062 TOULOUSE cedex 9

Dr Catherine FEUILLOLAY
Assay Manager

Dr Jocelyne BACARIA
Quality Manager

1. Test Laboratory

FONDEREPHAR

Faculté des Sciences Pharmaceutiques, 35 chemin des Maraîchers 31062 Toulouse cedex 9, France

2. Identification of the aerial disinfection system

Device : NOCOSPRAY

Serial number: 37S347

Disinfectant : NOCOLYSE®

Batch : 041017N (Expiry date: 10/2019)

Internal Code. 17-2217-1

Concentration of product in the room : 3 mL/m³

One treatment with one hour of wait (carriers recovery after waiting)

Amount of disinfectant diffusion \approx 96 mL/treatment (exactly 103,5g)

Promotor : OXY'PHARM

Storage conditions: Ambient temperature

Period of testing: November - December 2017

Actives Substances: Hydrogen peroxide (6%)

3. Experimental Conditions

a. Tests micro-organisms


- Sporicidal activity :
 - o *Clostridium difficile* (spore) NCTC 13366

The test surface is a blue curtain which was transferred by OXY'PHARM.

Reference : "Polypropylène incorporé à un agent antimicrobien à base d'argent"
Batch 036111 - AMG MEDICAL INC

b. Conditions of aerial disinfection system use

- Room :



Relative humidity ranging from 42% to 38%.

Initial temperatures ranging from 19,6°C to 19,2°C.

Test room volume : 32m³.

Distance between the appartus and the carriers : 2,6m (tableau B.1).

Squares (3 cm) have been drawn on the surface (rough side) at a height between 1,3m and 1,5m, in a vertical position (the curtain has been set up in the room), with facing away from the source

Diluants and culture media

Interfering substances

1/20 reconstituted milk (Internal preparation - Batch 6279 Exp. Dec/29/2017)

Diluants

Suspension preparation : EPPI (Cooper, Batch 19KL09GA Exp. October 2019)

Recovery solution (Internal preparation - Batch 7594 Exp. Dec/28/2017)

Filtration membranes

Nitrocellulos membranes 0,45µm (Millipore, batch F7JA35599 Exp. July 2019)

Culture media

Medium for *Clostridium difficile* (Internal preparation - Batch 7589 Exp. Dec/12/2017)

4. Assays

- **Treatment 3 mL / m³ - waiting 1H**

Tests microorganisms	N Test suspension (CFU/mL)	Preliminary assay			T Control (CFU/spot - 50µL)	n'1 + n'2 UFC/ spot 50µL (dilution/filtration - disc in agar)	Log reduction - Mean
		n1/N1	n2/N2	n3/N1			
<i>C. difficile</i>	2.10 ⁵ - 5.10 ⁵	n1 > 0.5 N1	n2 > 0.5 N2	n3 > 0.5 N1	≈ 10 ⁴		
Date Nov/30/2017	2,3.10 ⁵	d1 : 56/110 d2 : 55/110	d1 : 82/137 d2 : 89/137	d1 : 118/110 d2 : 101/110	d1 : 9,5.10 ² d2 : 1,6.10 ³	d1 : 0 + 0 d2 : 0 + 0	R1 : 3,11 R2 : 3,11 R = 3,11
B: 19,6°C/RH 42%							
E: 19,2°C/RH 38%					T = 1,3.10 ³ *		

* A 1 log decrease after the step of drying has been observed

T: counting of micro-organisms on the discs.

N₁ : counting of test suspension by pour plate technique - N₂ : counting of test suspension by filtration method

n₁ : counting to search inhibitor effect in agar medium - n₂ : counting to search inhibitor effect on membrane filtration - n₃ : counting to search inhibitor effect after inclusion of disc in agar medium

n'1 : number of survival micro-organisms in 100mL of tryptone-salt - n'2 : number of micro-organisms after inclusion of the carrier surface in agar medium.


n'1 + n'2 : total number of survival micro-organisms on the carrier surface.

d1 : disc N°1 / d2 : disc N°2 / d3 : disc N°3

FONDEREPHAR

Faculté des Sciences Pharmaceutiques - 35 Chemin des Maraîchers - 31062 TOULOUSE Cedex 09
Tél. 05 62 25 68 60 Fax. 05 61 25 95 72 Email. contact@fonderephar.com

No modification of blue color have been observed after the treatment



5. Conclusion

According to standard NF T 72-281 (November 2014) and for the test conditions, the couple NOCOSPRAY (serial number 37S347) / NOCOLYSE® (Batch 041017N Expiry date: December 2019), in use of human health, on the blue curtain transfer by client, led to:

- A sporicidal activity (log reduction ≥ 3)
 - After 3 mL/m³ treatment - 1 hour of wait on the following strain :
 - *Clostridium difficile* (spore) NCTC 13366