



Toulouse, July 7<sup>th</sup> 2023

## ASSAY REPORT N° 23-2081

### STUDY 23-3455

**STANDARD NF EN 17272 (Avril 2020)**  
**Chemical disinfectants and antiseptics -**  
**Methods of airborne room disinfection by automated process - Determination of**  
**bactericidal, mycobactericidal, sporicidal, fungicidal, yeasticidal,**  
**virucidal and phagocidal activities**

**Medical area**  
**Clean conditions**

**Additional microorganisms:** *S. aureus* ATCC 33591 (MRSA) and *P. aeruginosa* CIP 106817 (MDRP)

**Client**

**Company registration SIREN 448974253**  
**829 rue Marcel Paul**  
**94500 CHAMPIGNY SUR MARNE**  
**FRANCE**

**Assay laboratory**

**FONDEREPHAR**  
**Faculté des Sciences Pharmaceutiques**  
**35 Chemin des Maraîchers**  
**31062 TOULOUSE cedex 9**  
**FRANCE**

**Dr Catherine FEUILLOLAY**  
Study Manager

**Dr Jocelyne BACARIA**  
Quality Manager

## 1. Test Laboratory

FONDEREPHAR  
Faculté des Sciences Pharmaceutiques  
35 chemin des Maraîchers  
31062 Toulouse cedex 9  
France

## 2. Identification of the aerial disinfection system

Device : Diffuser PX-00  
Serial number : 172X731

Disinfectant : Formula N-1  
Batch : A310523N/1  
Exp.: May/2025  
Receipt : Jun/26/2023  
Internal code: 23-3455-1

Concentration of product: 8mL/m<sup>3</sup>

One treatment - Waiting time 120 minutes after the end of diffusion

Amount of disinfectant diffusion: 222g (difference of weight's flask between T<sub>0</sub> and at the end of the diffusion and the waiting time)

Time of diffusion : 15 minutes

Promotor : Company registration SIREN 448974253

Storage conditions: Ambient temperature

Period of testing: June - July 2023

Active Substance: Hydrogen Peroxide (6%)

**NB. Obligatory conditions of the standard (efficacy and distribution) in the test report N° 22-1907 (dated May/29/2022)**

## 3. Experimental Conditions

### a. Tests micro-organisms

- Bactericidal activity :
  - *Staphylococcus aureus* ATCC 33591 (MRSA)
  - *Pseudomonas aeruginosa* CIP 106817 (MDRP)


### b. Carriers

The selected tests surfaces are stainless steel discs, flats, corresponding to the requirements of paragraph 5.2.3.2 of the standard. The supplier is MERCIER CLAUSSE (France).

#### 4. Efficacy tests

##### a. Conditions of aerial disinfection system use

- Room :



Relative humidity ranging from 50% (see results).  
Initial temperatures ranging from 20,9°C (see results).

Test room volume : 32m<sup>3</sup>.

Distance between the appartus and the carriers : 2,60m (tableau B.1), 1,15m from floor.

##### b. Diluants, culture media and membranes

###### Interfering substances

1/20 reconstituted milk (Internal preparation - Batch 11660 Exp. Jul/26/2023)

BSA fraction V 0,3g/l (Internal preparation - Batches 577 and 579)

###### Diluants

Suspension preparation : Water for Injectable Preparations (WIP) (interference of product with Tryptone-salt) (Cooper - Batch 19QIAEW0 Exp. Aug/2024)

Recovery solution + 0,5% Tween 80 (Internal preparation - Batch 11601)

###### Filtration membranes

White nitrocellulose membranes 0,45 µm (Millipore - Batch F2DB43464)

###### Culture media

Trypcase soy agar (Biomérieux - Batch 1009871570 Exp. Aug/14/2024)

c. Results

- 8 mL / m<sup>3</sup> - waiting 120 minutes - Batch A310523N/1

Tests microorganisms	N Test suspension (CFU/mL)	Preliminary tests			T Control (CFU/spot - 50µL)	n'1 + n'2 CFU / spot 50µL (dilution/filtration - disc in agar)	Log reduction - Mean
		n1/N1	n2/N2	n3/N1			
	5.10 <sup>7</sup> - 2.10 <sup>9</sup>	n1 > 0.5 N1	n2 > 0.5 N2	n3 > 0.5 N1	≈ 10 <sup>6</sup>		
<i>S.aureus</i> ATCC 33591 Assay Jun/29/23 20,9°C/RH 50%	8,75.10 <sup>8</sup>	d1 : 96/88 d2 : 91/88	d1 : 91/84 d2 : 78/84	d1 : 82/88 d2 : 78/88	d1 : 2,04.10 <sup>7</sup> d2 : 2,12.10 <sup>7</sup>  T = 2,08.10 <sup>7</sup>	E1 : 0 + 0 E2 : 0 + 0 E3 : 1 + 0	R1 : 7,32 R2 : 7,32 R3 : 7,32 <b>R = 7,32</b>

T: counting of micro-organisms on the discs.

N<sub>1</sub> : counting of test suspension by pour plate technique - N<sub>2</sub> : counting of test suspension by filtration method

n<sub>1</sub> : counting to search inhibitor effect in agar medium - n<sub>2</sub> : counting to search inhibitor effect on membrane filtration - n<sub>3</sub> : counting to search inhibitor effect after inclusion of disc in agar medium

n'1 : number of survival micro-organisms in 100mL of tryptone-salt - n'2 : number of micro-organisms after inclusion of the disc in agar medium.

n'1 + n'2 : total number of survival micro-organisms on the carrier surface.

d1 : disc N°1 / d2 : disc N°2 / d3 : disc N°3

- 8 mL / m<sup>3</sup> - waiting 120 minutes - Batch A310523N/1

Tests microorganisms	N	Preliminary tests			T	n'1 + n'2 CFU / spot 50µL (dilution/filtration - disc in agar)	Log reduction - Mean
	Test suspension (CFU/mL)	n1/N1	n2/N2	n3/N1	Control (CFU/spot - 50µL)		
	5.10 <sup>7</sup> - 5.10 <sup>9</sup>	n1 > 0.5 N1	n2 > 0.5 N2	n3 > 0.5 N1	≈ 10 <sup>6</sup>		
<i>P.aeruginosa</i> CIP106817 Assay Jun/29/23 20,9°C/RH 50%	3,8.10 <sup>9</sup>	d1 : 27/38 d2 : 31/38	d1 : 38/31 d2 : 34/31	d1 : 30/38 d2 : 28/38	d1 : 7,5.10 <sup>6</sup> d2 : 2,84.10 <sup>6</sup>  T = 5,17.10 <sup>6</sup>	E1 : 2 + 0 E2 : 67 + 0 E3 : 5 + 0	R1 : 6,41 R2 : 4,88 R3 : 6,01 <b>R = 5,77</b>

T: counting of micro-organisms on the discs.

N<sub>1</sub> : counting of test suspension by pour plate technique - N<sub>2</sub> : counting of test suspension by filtration method

n<sub>1</sub> : counting to search inhibitor effect in agar medium - n<sub>2</sub> : counting to search inhibitor effect on membrane filtration - n<sub>3</sub> : counting to search inhibitor effect after inclusion of disc in agar medium

n'1 : number of survival micro-organisms in 100mL of tryptone-salt - n'2 : number of micro-organisms after inclusion of the disc in agar medium.

n'1 + n'2 : total number of survival micro-organisms on the carrier surface. d1 : disc N°1 / d2 : disc N°2 / d3 : disc N°3

## 5. Conclusion

The device/product combination: diffuser PX-00 serial number 172X731 / Formula N-1 (batch A310523N/1), for use in clean conditions, in medical area, meets the criteria of standard NF EN 17272 (April 2020) for bactericidal efficacy test with additional microorganisms (*S. aureus* ATCC 33591 (MRSA) and *P. aeruginosa* CIP106817 (MDRP)) after treatment at 8 mL/m<sup>3</sup> and a waiting time of 120 minutes.

The results hold only for the device/product under assay and apply to the sample as received.