



# The Oxy'Pharm biodegradable solution

An efficient and eco-responsible approach to combat all infections and contamination in scientific animal houses.

Hygiene is a major issue in scientific animal houses. Oxy'Pharm devices meet your needs in terms of versatility and guarantee impeccable bio-cleaning and bio-disinfection.

**Because the solution is always a combined one**, Oxy'Pharm has implemented a complete hygiene concept with its range of Sanivap & Nocotech products for external disinfection. The Alpha Clean range provides bio-cleaning of your instruments before the sterilisation process.

A perfectly adapted solution for scientific animal houses (sterilisation areas, airlocks, housing and food storage areas, etc.) because since the treatments are biodegradable, they have no impact on electronic equipment and are harmless to health.

₩ We offer all of our customers training in the use of our Sanivap and Nocotech devices.

BIO-CLEANING OF SURFACES

BIO-DISINFECTION OF SURFACES

BIO-CLEANING OF INSTRUMENTS







#### FOR OPTIMAL HYGIENE IN ALL YOUR SPACES...



















## **Steam bio-cleaning**

Leader in steam bio-cleaning, the Sanivap range has steam generators (SP400, SP540H, SP500/600 and SP800) and a variety of treatment accessories for thorough and effective cleaning.

Sanivap's devices ensure the bio-cleaning of surfaces on a daily basis: premises, furniture, equipment, etc.

#### **ADVANTAGES**



#### **EFFICIENCY AND SAFETY**

Effective on multi-resistant micro-organisms and safe for human use.



#### **ECONOMICAL**

The steam method no longer requires the use of detergents and disinfectants and reduces water consumption.



#### **STEAM QUALITY**

High-temperature steam, under pressure and with disinfectant while constantly controlling the parameters.



#### **RELIABILITY AND PERFORMANCE**

Intensive professional use and battery life of 8 to 10 hours per day.



#### **EASE OF USE**

Mobile cart, simple user interface, water level and control indicator.

## PRODUCTS ADAPTED TO YOUR NEEDS



#### **THE SP400 KIT**

Our 2850 W steam generator is specially designed to be used up to 8 hours a day.

Integrated into its bio-cleaning trolley, its use is simple and particularly well-suited to the daily maintenance protocols of animal houses.

Sanivap devices are Class IIa medical devices - Notified body: GMED 0459 - Manufacturer: Sanivap France

#### **ACCESSORIES**















## 100% automated bio-disinfection

Nocotech is an automated concept for surface disinfection by dry diffusion. The range is based on the combined action of an automated surface treatment device (Nocospray) and a liquid disinfection product (Nocolyse range). Insect control treatment can be done with the Oxypy range. Aerial diffusion can be used for the treatment of all accessible surfaces.

#### **ADVANTAGES**



#### **NO RISK**

Biodegradable, residue-free, non-corrosive (including on electronic surfaces), non-toxic, non-allergenic.



#### **ECONOMICAL**

Very low quantity of product required.



#### **EFFICIENT AND CERTIFIED**

Bactericide, virucide, sporicide and fungicide (validated by standard EN 17272).



#### **INSECTICIDE**

Insecticide treatment against scabies, bedbugs and crawling and flying insects.




#### **GUARANTEED TIME SAVING**

No action required before treatment (except cleaning) or after treatment (no need to rinse).

## PRODUCTS SUITED TO YOUR NEEDS

#### **NOCOSPRAY**

Our lightweight, portable surface biodisinfection machine. Volume of the room to be treated from 10 to 1000 m<sup>3</sup>.





### WALL MOUNTED NOCOSPRAY

Our surface bio-disinfection machine is particularly suitable for airlocks.

Volume of the room to be treated from 10 to 1000 m<sup>3</sup>.



#### **NOCOLYSE FOOD**

Our silver-free ready-to-use disinfectant with 7.9% hydrogen peroxide.

Available in 1L bottle, 5L can, 20L can.



 $\label{thm:condition} \textbf{Use biocidal products carefully. Read the label and product information before use.}$ 

#### **ACCESSORIES**



#### **REMOTE CONTROL**

For activation of the Nocospray from outside the room.



#### **PROGRAMMER**

For delayed activation of the Nocospray.



#### WALL CROSSING SYSTEM KIT

For disinfecting a room without entering it.



#### **DOUBLE NOZZLE**

For two-way spraying.



# Instruments pre-disinfection

The Alpha Clean range consists of an auto-washer and a neutral product (Clean One solution) which provides the automatic and complete treatment of all instruments after use and before sterilisation. The device has a single tray that opens like a drawer. Soiled instruments from the treatment area, are placed in there directly. The duration of the cycle is 45 min.

#### **ADVANTAGES**



#### **NO RISK**

With soiled instruments no longer requiring handling, any risk of infection for healthcare staff is eliminated.



#### **AUTOMATIC DOSING**

The disinfection/cleaning product is automatically dosed during each cycle by a peristaltic pump.



#### **EFFICIENT AND CERTIFIED**

The ultrasonic system disinfects the hardest to reach corners of the most complex instruments.



#### **ECONOMICAL**

Very small quantity of product required.



#### **GUARANTEED TIME SAVING**

No handling required before or after treatment.

## PRODUCTS ADAPTED TO YOUR NEEDS

#### **ALPHA CLEAN**

Our instrument disinfection device with ultrasonic tank.

Available in 16 L and 27 L capacities.



#### **CLEAN ONE**

Our pre-disinfectant liquid detergent for temperature stable and temperature sensitive instruments. It is biodegradable, bactericidal, yeasticidal and virucidal.



Available in a 5 L container.

#### **HOW DOES IT WORK?**

#### 1. PRE-DISINFECTION PHASE

The instruments placed in the tank are immersed in the Clean One disinfection solution.

#### 2. CLEANING PHASE

Ultrasound in combination with the Clean One solution disinfects every last corner of the instruments.

#### 3. RINSING PHASE

The instruments are immersed in clear water and the ultrasound action provides efficient rinsing.

#### **4. DRYING PHASE**

The turbine generates hot air up to 80°C.

#### **ACCESSORIES**

#### BASKET FOR ALPHA CLEAN 27 L

Provides you with a spare basket.

#### BASKET FOR ALPHA CLEAN 16 L

Provides you with a spare basket.

#### **DRILL BASKET**

Provides the option for several drills to be placed in the treatment device.

#### **CANNULA STAND**

Allows an optimised layout of the suction cannulas.

### STACKABLE BASKETS

Provides the option for several rows of instruments.

## **Notes**



www.oxypharm.net @ 829 rue Marcel Paul, 94500 Champigny-sur-Marne, France  $\boxtimes$  info@oxypharm.net % +33 1 45 18 78 70